

# Low-Cost Determination of Anions in Municipal Drinking Water Using EPA Method 300.0 and the ICS-900 IC System

## **INTRODUCTION**

U.S. EPA 300.0<sup>1</sup> is an ion chromatography (IC) method approved by the United States Environmental Protection Agency (EPA) to determine inorganic anions (fluoride, chloride, nitrite, sulfate, bromide, nitrate, and phosphate) in municipal drinking water and wastewater. Three of these anions (fluoride, nitrite, and nitrate) are considered contaminants in drinking water under the EPA's National Primary Drinking Water Regulations. High levels of fluoride cause bone disease; nitrite and nitrate can cause birth defects. Chloride and sulfate are not harmful to public health but should not exceed 250 mg/L for the proper taste of water, according to the EPA's National Secondary Drinking Water Regulations. The EPA enforces these regulations to ensure use of approved analytical methods when analyzing water samples to meet federal monitoring requirements or to comply with drinking water regulations.

The work shown here describes the determination of inorganic anions in drinking water with the ICS-900 and EPA method 300.0. The ICS-900 is an integrated, single-channel IC system designed for isocratic applications with suppressed conductivity detection.

The system facilitates continuous chemical suppression with MMS™ MicroMembrane™ Suppressor technology to enhance detector response, provide low baseline noise, and prevent baseline drift. The ICS-900 is designed for low maintenance costs and maximum uptime, making it a great choice for routine analysis.

To demonstrate the capability of the ICS-900 for water analysis, municipal drinking water of the city of San Jose, CA was analyzed at ambient temperature with the IonPac® AS22 column. The experimental result is shown in Figure 1.

All seven anions were separated and eluted within 13 min. Chloride was the most abundant anion in the sample, followed by sulfate, and their concentrations were determined to be 98 and 48 mg/L, respectively. Fluoride and nitrite were detected at less than 1 mg/L, which was below the maximum amount allowed by EPA regulations. The concentration of nitrate was determined to be 2.43 mg/L, which met the safety criterion (less than 10 mg/L) set by the National Primary Drinking Water Regulations.

Now sold under the  
Thermo Scientific brand

**Thermo**  
SCIENTIFIC

## CONDITIONS

Dionex ICS-900, Dionex AS-DV Autosampler, and Chromeleon® Chromatography Data System (CDS) software are used in this experiment. All experimental parameters are listed in Figure 1.

## SAMPLE INFORMATION AND PREPARATION

Municipal drinking water of the city of San Jose, CA was analyzed by IC without sample pretreatment.

## CONCLUSION

The Dionex ICS-900 provides a complete solution for routine analysis of inorganic anions in water and meets all requirements specified in EPA method 300.0. A thorough discussion on anion determinations in municipal drinking water using carbonate eluents can be found in AN 133.<sup>2</sup>

## REFERENCES

1. *The Determination of Inorganic Anions in Water by Ion Chromatography; Method 300.0*, Revision 2.1; U.S. Environmental Protection Agency; Cincinnati, OH, 1993.
2. Dionex Corporation, *Determination of Inorganic Anions in Drinking Water by Ion Chromatography*. Application Note 133, LPN 1192, Sunnyvale, CA, 2004.




Figure 1. Separation of anions in municipal drinking water sample on the IonPac AS22 column using the ICS-900 system.

IonPac and Chromeleon are registered trademarks, and MMS and MicroMembrane are trademarks of Dionex Corporation.

Passion. Power. Productivity.



### Dionex Corporation

1228 Titan Way  
P.O. Box 3603  
Sunnyvale, CA  
94088-3603  
(408) 737-0700

### North America

U.S./Canada (847) 295-7500

### South America

Brazil (55) 11 3731 5140

### Europe

Austria (43) 1 616 51 25 Benelux (31) 20 683 9768 (32) 3 353 4294  
Denmark (45) 36 36 90 90 France (33) 1 39 30 01 10 Germany (49) 6126 991 0  
Ireland (353) 1 644 0064 Italy (39) 02 51 62 1267 Sweden (46) 8 473 3380  
Switzerland (41) 62 205 9966 United Kingdom (44) 1276 691722

### Asia Pacific

Australia (61) 2 9420 5233 China (852) 2428 3282 India (91) 22 2764 2735  
Japan (81) 6 6885 1213 Korea (82) 2 2653 2580 Singapore (65) 6289 1190  
Taiwan (886) 2 8751 6655

[www.dionex.com](http://www.dionex.com)



LPN 2711 PDF 1/11  
©2011 Dionex Corporation