




# Sample Preparation for Ion Chromatography

Comprehensive Solutions to Sample Preparation Challenges



Now sold under the  
Thermo Scientific brand

**Thermo**  
SCIENTIFIC

 **DIONEX**

Speed · Simplicity · Solutions

# Choose from a Broad Range of Solutions


No matter what your needs are, turn to Dionex—the ion chromatography (IC) technology leader—for a comprehensive range of proven sample preparation solutions.

- Easy to understand and implement
- Simple to operate
- Highly reliable and reproducible
- Saves time and lowers operating costs


The simplicity and ruggedness of IC often requires no need to prepare a sample before injection. In many cases, a simple dilution of the sample with deionized water is all that's needed. Sometimes more is required to protect the column by removing contaminants or particles or to enhance the sample to make it detectable. Any sample preparation technique must maintain sample integrity, be reproducible, and require minimal labor to accomplish.

## Dionex solutions span:

- Off-line solutions that are easy and practical to implement
- In-line and automated technologies using proprietary technologies
- Specialized and high-throughput solutions



Autosampler Solutions



In-line



RFIC-ESP



Chromeleon Software



Cartridge Technologies



Off-line



Electrolytic  
Trap Devices



IC Integration

# Dionex Sample Preparation Solutions

## AutoDilution

Conditional determination of high analyte concentration with automatic dilution and reinjection

- Autosampler with Chromeleon® Chromatography Data System (CDS) Software

## Matrix Elimination

Remove interfering matrix ions and organic contaminants

- Off-line OnGuard® II Cartridges
- In-line InGuard® Cartridges
- In-line IonPac® NG-1 Columns

## Concentrator Columns

Capture and concentrate ions from large, dilute sample volumes

- 15 Anion IC and RFIC Concentrator Columns
- 3 Cation IC and RFIC Columns
- 3 Transition Metal IC and RFIC Trap Columns

## Trap Columns

Purify eluent or polish samples

- Anion IC and RFIC Trap Columns
- Cation IC and RFIC Concentrator Columns

## RFIC-ESP

Reagent-Free™ Ion Chromatography (RFIC™)-Electrolytic Sample Preparation (ESP)

- Continuously Regenerated Trap Columns (CR-TC)
- AutoNeutralization™ Cartridges
- Electrolytic Water Purifiers (EWPs) for Trace Analysis
- AutoPrep with ESP for Trace Analysis of Anions or Cations

## Sample Filtration

Remove particulates that can damage your system; remove microbes to prevent analyte degradation

- Off-line Syringe Filters
- In-line AS-DV Autosampler Filters
- In-line Filter Frit Cartridges

## Specialty Systems

Combustion IC and Air Sampling systems

- Combustion IC for Halogen and Sulfur Analysis
- Air Sampling System for Ionic Analysis

## Chromeleon Chromatography Data System Software

Plays an integral role in most Dionex automated solutions

- Monitor system pressure to detect filter clogs
- Optimum injection parameters for injection loops and concentrator columns
- Visualization tools for matrix elimination and analyte concentration
- Conditionally controlled autodilution and electrolytic sample preparation techniques

# AutoDilution and Analysis

## Automatic Detection, Dilution, and Re-injection of Highly Concentrated Samples

### Challenge

- Sample too concentrated for analysis
- Sample needs dilution
- Delay in analysis time
- Extra work

### Solution—Dionex AutoDilution Capabilities




- Automated detection of concentrated samples
- Automatic dilution, reinjection, and analysis
- Faster results
- Saves effort, time, and money

IC systems equipped with an autosampler, Chromeleon CDS software, and the AutoDilution license can monitor chromatographic results. If peak results fall



outside a user-specified range, the system will automatically dilute the sample and reinject it. Decisions to dilute and reinject can be based on a single parameter such as

peak area, height, amount, resolution or other chromatographic value, or on a combination of values.


### AutoDilution: Large Loop/Small Loop Conditional Injection



Sample is too concentrated, and analytes are out of the calibration range.



Data system detects concentrated sample, and then autosampler dilutes and reinjects it. Analytes are now in the calibration range.



# Matrix Elimination

## Selectively Remove Interfering Ions from Matrices Using Solid-Phase Extractants

### Challenge

Interfering ions  
 Poor IC performance  
 Short consumables lifetime  
 High cost of ownership

### Solution—Dionex Matrix Elimination Cartridges



Eliminate interfering ions  
 Automated and manual methods  
 Superior IC performance  
 Extended lifetime of consumables  
 Lower cost of ownership

### Off-line Matrix Ion Elimination


#### OnGuard II Cartridges

OnGuard II Cartridges have Luer inlets for manual processing of a limited number of samples. A wide range of cartridges loaded with resins in 1 and 2.5 cc capacity formats remove matrix interferences encountered in many ion chromatography applications, (see table at right).

These cartridges have wide pH stability and allow low-level ion analysis. They remove both column-fouling substances and interfering ions where several mL of sample is available. An OnGuard Sample Prep Station is available for simultaneous matrix elimination of multiple samples.

OnGuard II Cartridges*	
Resin	Removes from Sample Matrices
A	Anionic contaminants; Neutralization of highly acidic samples
Ag	Chloride, bromide, and iodide
Ba	High concentrations of sulfate
H	High levels of alkaline earth and transition metals; Neutralizes highly alkaline samples (e.g., sodium hydroxide or sodium carbonate)
Na	High levels of alkaline earth and transition metals without acidifying the sample
M	Transition metals; also for matrix elimination of alkali and alkaline earth metals
P	Phenolic fraction of humic acids, tannic acids, lignins, anthocyanins, and azodyes prior to analysis by ion exchange
RP	Hydrophobic substances such as aromatic dyes, some aromatic carboxylic acids, hydrocarbons, and surfactants
Ag/H	Chloride, bromide, and iodide from concentrated matrices such as brines
Ba/Ag/H	Respective contaminants using a layered resin cartridge containing the Ba, Ag, and H styrene-based sulfonic acid resins

\*Cartridges can be used singly or in series.



OnGuard II RP for the removal of fat from milk. Failure to remove fat prior to IC analysis leads to high column backpressures, loss in column capacity, and decreased retention times.

# Matrix Elimination


## Continued

### In-line Matrix Ion Elimination InGuard Cartridges

InGuard Cartridges provide high-throughput automated removal of matrix interferences encountered in many IC applications. Depending on the cartridge's location in the

IC system, the autosampler (low pressure) or pump (high pressure) pushes the sample through the InGuard Cartridge. The number of samples that can be treated with an InGuard Cartridge depends on the matrix composition, concentration of compounds to be removed,

and the sample injection volume. Chemistries operate optimally at flow rates between 0.5 and 2 mL/min. They are suitable for both standard-bore and microbore IC systems.



Treatment of a brine sample with InGuard Ag Cartridge removes interfering chloride ions to permit analysis of small amounts of other anions.

### InGuard Cartridges for High-Throughput Automated Matrix Elimination\*

Resin	Removes from Sample Matrices
Ag	Halides
H	Cations and pH adjustment
Na	Cations
HRP	Hydrophobic species, azo- and cyano-containing species
Na/HRP	Metals (Na) plus organics (HRP)

\*InGuard Kits are available that contain everything you need for high-throughput matrix elimination of organics, halides, sample pH reduction, and more.

### IonPac NG-1 Columns

These columns are packed with highly hydrophobic cross-linked polystyrene resin beads for the removal of organic

contaminants from samples, such as humic acids. The option of a two-column setup with a 10-port valve enables automatic resin cleaning of one trap column while the other

is being used for sample preparation. NG-1 kits are available for organic removal for anion and cation analysis.


# Concentrator Columns

## Capture and Concentrate Sample Ions from Large Sample Volumes

### Challenge


- Sample ions too dilute
- Manual methods prone to errors and not reproducible
- Time consuming

### Solution—Dionex Sample Pre-Concentration



- Preconcentrate ions, increase sensitivity
- Fully automated, highly reproducible
- Reduced errors and contamination
- Saves time and money

IonPac Concentrator Columns can be used to preconcentrate ions as required for high-purity water analysis in, for example, the power and semiconductor industries. The columns retain ions from a measured volume of aqueous sample matrix, thereby concentrating the analyte species and lowering detection limits. Concentrator Columns enable the analysis of ions at  $\mu\text{g/L}$  (ppb) to  $\text{ng/L}$  (ppt) levels.



The concentrator column is typically placed on the injection valve in place of the sample loop, just before the analytical column. Sample is first collected on the concentrator column, which strips and retains the ions as the sample flows to waste. With a switch of the valve, eluent then flows through the column eluting the analytes to the analytical column.

### Anion IC and RFIC Concentrator Columns


Dionex offers 15 Anion Concentrator Columns with differing sizes and capacities for applications using hydroxide, carbonate/bicarbonate, or borate eluents.

### Cation IC and RFIC Concentrator Columns

Cation Concentrator Columns are fully functionalized with carboxylic acid for methane sulfonic acid (MSA) and sulfuric acid eluents. One column has a sulfonic acid functional group for  $\text{HCl/DAP}\cdot\text{HCl}$  eluents.

### Transition Metal IC and RFIC Concentrator Columns

The Trace Metal Concentrator Column is a high-capacity, cation concentrator column used for coupling the MetPac CC-1 column to the IonPac CS5 or CS5A analytical columns when performing chelation IC.



Comparison of direct injection and pre-concentration using the UTAC-LP1 for the determination of trace anions.

# Trap Columns

## Purify Eluent or Polish Samples

### Challenge

Unwanted analytes interfering with separation

Contaminant buildup on guard and/or analytical column

Poor results

### Solution—Dionex Trap Columns




Remove eluent and sample contamination and improve results

Easy to use

Prolong lifetime of consumables

Lower cost of ownership

IonPac Trap Columns prevent unwanted analytes from interfering with the separation of your analytes of interest. They contain high-capacity, low-efficiency, ion-exchange resins that strip and retain trace contaminants from the eluent and prevent their concentration on the guard and analytical columns.



CR-CTC II neutralizes acidic samples, removes anion contaminants.

### Trap Columns and RFIC Devices

These columns contain high-capacity resins designed to efficiently trap contaminants in eluents or samples. Back pressure is strictly controlled to provide the most useful devices available. Included are ATC-3 and ATC-HC anion trap devices, as well as the borate form ATC column. The CTC Cation Trap column provides high efficiency trapping of cations from samples, and in the hydrogen form it acts to neutralize basic samples.

### Continuously Regenerated Trap Columns (CR-TC)

CR-TC technology provides continuously regenerated sample purification devices for maintenance-free operation. Both anion removal (CR-ATC) and cation removal (CR-CTC II) devices are available. CR-CTC II removes lithium and neutralizes hydroxide of borated water samples from nuclear power plants.

### MFC-1 Metal-Free Trap Columns

The IonPac MFC-1 Metal-Free Column contains a special chelating resin. This column is placed in the eluent line prior to the injection valve to remove transition metals from high-pH eluents.

### Cation Polisher Columns for Anion Analysis

These columns remove metallic contaminants and other cations such as calcium and magnesium from the sample stream while performing anion analysis, and include the:

- Cation Polisher  $\text{Na}^+$  Column for anion analysis, which is a cation-exchange column in the sodium form and packed in a low-pressure format for autosampler operation.
- Cation Polisher  $\text{H}^+$  Column for anion analysis, which is a higher-capacity version in the hydronium form.



# RFIC-ESP

## Dionex Reagent-Free Ion Chromatography Electrolytic Sample Preparation

### Challenge

- Ionic contamination in eluent stream
- Matrix contamination
- Cross-contamination for trace analysis
- Method difficult to transfer between systems

### Solution—RFIC-ESP Options



- Automated ion removal from eluent stream
- Easy to understand
- Simple transfer to another system
- Highly reproducible and precise
- Reduced handling and human error

### Electrolytic Devices for Automated Sample Pretreatment

Several Dionex IC systems integrate hardware and software to enable automated electrolytic removal of ions from an eluent stream using protocols that are easily understandable, simple to transfer to another system, and highly reproducible and precise.

- Highly reproducible results, instrument-to-instrument and lab-to-lab
- Reduced human error and inconsistencies
- Reduced manual operation with reduction in labor costs and sample processing time
- Reduced system complexity by eliminating sample loading pump
- Ultralow backgrounds from on-line electrolytic water purification with routine trace (ng/L) analysis
- Automated purification, concentration, and calibrations providing reduced labor costs and improved quality of results

### AutoNeutralization Cartridges

The determination of anions in concentrated bases and the determination of cations in concentrated acids are important applications in industrial, health and safety, and environmental fields. Direct injection of concentrated acids or bases will overload

the column resulting in poor chromatography and making ion quantification difficult or impossible. AutoNeutralization is a high-capacity electrolytic ion-exchange device that neutralizes basic or acidic samples prior to injection. A sample can be trapped in the AutoNeutralization device for as long as necessary to fully neutralize the sample matrix. The sample is then transferred to a concentrator column, and then to the analytical separation column.

### Electrolytic Water Purifiers (EWPs)

EWP units are devices packed with various ion-exchange resins, and produce low flow rate streams for water of exceptional purity. These water streams can be used to displace water for loading samples onto concentrator columns, either directly or via another sample preparation device such as the CR-TC.

### RFIC-ESP and AutoPrep

AutoPrep on the ICS-2100 Ion Chromatography System utilizes RFIC-ESP techniques that simplify the analysis of trace levels of anions or cations. The use of EWPs to transfer large sample volumes to the concentrator column virtually eliminates background contamination, allowing for low ng/L concentrations of ions to be conveniently

analyzed. Proportional volume loops enable injection of standards of concentrations 1000 times higher than samples. Automated multilevel calibration is easily accomplished using Chromeleon software. Anion and cation determinations as low as 1 ng/L (ppt) are easily achievable.

### RFIC-ESP Advantageous for:

- Trace analysis of ions in power plant and semiconductor ultrapure water samples
- Determination of ions in concentrated acids and bases after neutralization
- Cation elimination from power plant borated waters with continuous regeneration of the trap device
- In-line neutralization of dilute samples
- Continuous desalting of an ion-exchange eluent for coupling to mass spectrometers
- Continuous desalting of ion-exchange eluents before fraction collection


# Sample Filtration

## Remove Particulates to Protect Valuable Chromatography System Components

### Challenge

- Particulate blockage in system
- High backpressure
- Instrument downtime
- Microbial contamination
- Sample degradation

### Solution—Dionex Filtration Devices



- Eliminates column contamination and high system back pressures
- Manual and automated methods
- Reduces cost of ownership
- Improves reliability and reproducibility of results

### Off-line

#### Particle Filters

- Practical and easy to use
- Manually process samples with a 0.45 µm nylon membrane

The off-line particle filter has a sterile syringe Luer fitting for injecting a sample through a 0.45 µm nylon filter. This device is perfect for filtering a limited number of samples, such as collecting samples in the field to comply with EPA requirements to filter microbes.

### In-line

#### AS-DV Autosampler Vial Filters

- Automated, robust, inexpensive, and easy to use
- Efficient filtration for highly contaminated samples
- No sample carryover

Used with the AS-DV high-performance autosampler, these 20 µm filtering devices sit atop the sample in either the 5 mL or 0.5 mL sample vial for automatic, high-capacity removal of particulates prior to injection. Because particulates settle to the bottom of the vial, the least contaminated

sample is filtered first, thus preventing clogging and enhancing the filter's capacity. Individual filters are used for each sample, eliminating carryover concerns.

#### In-line Filters

- Housing made of PEEK™ prevents ion contamination
- Minimal dead volume (~1.3 µL) for use in various in-line locations
- High pressure rating of 5000 psi
- Several in-line filtering devices to match your application

Two types of in-line sample filters are available:

#### High-Capacity In-line Filter Cartridges

Used in the sample and/or eluent line for high-capacity filtration, this in-line filter cartridge has 35, 5, and 0.45 µm filters.

#### Low-Dead Volume In-line Filter Cartridges

These in-line filter cartridges have 0.5 µm PEEK filter frits and provide low dead-volume sample filtration and added system protection when low particulate matter is expected.

### Kits Add Convenience for In-line Filtering

- A variety of in-line filter kits are available with everything needed to configure your IC system.

# Specialized Systems

## Combustion Ion Chromatography

### Automated analysis of halogens and sulfur in petrochemicals, industrial materials, and consumer products

Combustion IC enables anion analysis in a variety of difficult-to-process samples such as ores, printed circuit boards, solders, plastics, and petroleum products, among others. The sample combustion process automates

conversion of a solid, liquid, or gaseous sample to an aqueous solution that can be injected directly into an IC system for fully automated sample processing. Samples are efficiently oxidized by heating to 1000 °C in an

air or oxygen stream for just a few minutes. Volatile combustion products are carried by the combustion gas stream into a sparger where they are trapped in water. The sample is then loaded onto the IC injection valve and analyzed.



The Combustion Ion Chromatography System comprises (from right to left) an autosampler, combustion unit, absorption unit, and ion chromatography system.

## Air Sampling System

### Automated analysis of volatile ionic constituents in air from emissions contributing to acid rain

Volatile ionic constituents in air result from emissions ultimately become acid rain. Particulate-borne ions can be dissolved in precipitation. The URG air sampling system easily conducts air monitoring on a

continuous basis. The system separates and collects particulate material and volatiles, and captures the volatile constituents using an aqueous sparging system. Particulates, dissolved in water automatically, and

samples are then automatically injected and analyzed by ion chromatography. The system can be operated for days without user intervention.



The URG Air Sampling System with integrated ion chromatographic analysis provide near real-time results.

## Enjoy Industry-Leading Support

Customer Support Centers for Dionex products are located in the Americas, Europe, and Asia. State-of-the-art laboratories are equipped with the full line of Dionex consumables, instruments, and software capabilities. Support Centers provide accessible locations for advanced training and enhanced application development capabilities. Users can attend these laboratories to learn new skills in addressing challenging applications, receive training and support, and discover new, innovative IC and sample preparation solutions.

©2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

<b>U.S./Canada</b>	(847) 295 7500	<b>Denmark</b>	(45) 36 36 90 90	<b>Sweden</b>	(46) 8 473 3380	<b>India</b>	(91) 22 2764 2735
<b>Brazil</b>	(55) 11 3731 5140	<b>France</b>	(33) 1 39 30 01 10	<b>Switzerland</b>	(41) 62 205 9966	<b>Japan</b>	(81) 6 6885 1213
<b>Austria</b>	(43) 1 616 51 25	<b>Germany</b>	(49) 6126 991 0	<b>United Kingdom</b>	(44) 1276 691722	<b>Korea</b>	(82) 2 2653 2580
<b>Benelux</b>	(31) 20 683 9768	<b>Ireland</b>	(353) 1 644 0064	<b>Australia</b>	(61) 2 9420 5233	<b>Singapore</b>	(65) 6289 1190
	(32) 3 353 42 94	<b>Italy</b>	(39) 02 51 62 1267	<b>China</b>	(852) 2428 3282	<b>Taiwan</b>	(886) 2 8751 6655