

# Determination of Inorganic Anions Using the Dionex IonSwift MAX-200 Monolith Column with Capillary Ion Chromatography

Fei Pang, Terri Christison, Barbara Shao, and Cathy Tanner  
Thermo Fisher Scientific, Sunnyvale, CA, USA

## Overview

**Purpose:** To demonstrate the fast determination of bromate in a municipal drinking water using the monolith-based Thermo Scientific Dionex IonSwift™ MAX-200 capillary column.

**Methods:** Bromate is separated from other inorganic anions using an hydroxide gradient and detected by suppressed conductivity detection.

**Results:** A minimum of 1 µg/L bromate can be detected. Good linearity and recovery of bromate are observed.

## Introduction

Ion chromatography (IC) with suppressed conductivity detection is a well-established technique for the determination of inorganic and organic ions in the water analysis industry. The recent introduction of high-pressure capillary IC brings additional advantages to the analysis of inorganic ions such as speed and high resolution. Due to low consumption of eluent, capillary IC systems can be operated continuously and are always ready for analysis, thus redefining the workflow for IC and improving method performance. The operation of capillary IC at higher pressures and/or higher flow rates improves separation efficiency and/or speed. This work shows the determination of inorganic anions and oxyhalides using the Dionex IonSwift MAX-200 capillary column.

## Methods

### Sample Preparation

Bromate at 10 µg/L was spiked into a municipal drinking water sample.

## Ion Chromatography Equipment and Data Analysis

Thermo Scientific Dionex ICS-5000 RFIC™ Capillary Ion Chromatography System consisting of:

- DP Dual Isocratic Capillary Pump
- DC Detector/Chromatography Module
- Thermo Scientific Dionex IC Cube™ Capillary Module (Figure 1)
- ICS-5000 CD Capillary Conductivity Detector for anions and cations
- EG Eluent Generator Module
- Thermo Scientific Dionex AS-AP Autosampler with diverter valve
- Thermo Scientific Dionex Chromeleon™ Chromatography Data System software

The Dionex ICS-5000 IC Cube contains the consumables with easy-to-install cartridges: the EG Degas module, column tray, four-port injection valve with 0.4 µL internal loop, suppressor module, and Carbonate Removal Device (CRD) module.


FIGURE 1. Dionex ICS-5000 IC Cube Module.



## Conditions

All the experimental conditions are listed Figures 4–8. A six-port injection valve with a 3  $\mu\text{L}$  sample loop (Figure 2) was used for all separations replacing the standard four-port injection valve and 0.4  $\mu\text{L}$  sample loop.

FIGURE 2. A six-port injection valve with 3  $\mu\text{L}$  sample.




## Results

### Inorganic Anions Separation on Dionex IonSwift MAX-200 Capillary Column

Dionex IonSwift MAX-200 capillary column is designed to provide high-speed, high-resolution inorganic anions separations using a hydroxide gradient delivered by an eluent generator. Speed and resolution are two competing performance factors in conventional chromatography separations using porous media. A solution for this problem is a new generation of separation media: monolithic structures. Monoliths contain a network of large separation pores where the mass transfer is achieved primarily by convective flow instead of diffusion during separation. The highly improved mass transfer at high-flow velocity in monoliths removes the peak broadening problem in conventional porous separation media. High-speed separations can be achieved with minimal loss of resolution. Therefore, monoliths offer high-resolution and high-speed separations similar to well constructed nonporous media. Figure 3 shows the structure of the Dionex IonSwift MAX-200 capillary column. Determination of inorganic anions at different flow rates is shown in Figure 4. Figures 5 and 6 demonstrate that inorganic anions can be determined within 10 min at 25  $\mu\text{L}/\text{min}$  flow rate taking advantage of the high-pressure capabilities of capillary IC and the relatively low back pressure provided by a monolith column.

FIGURE 3. Structure of Dionex IonSwift MAX-200 column.




29077

FIGURE 4. Fast determination of inorganic anions using the Dionex IonSwift MAX-200 capillary column.

Column: Dionex IonSwift MAX-200,  
0.25 mm  $\times$  250 mm  
Eluent Source: Dionex capillary EGC-KOH cartridge  
Gradient: KOH gradient from 2 mM to 50 mM  
Flow Rate: A: 10, B: 15, C: 20, D: 25  $\mu\text{L}/\text{min}$   
Inj. Volume: 0.4  $\mu\text{L}$   
Column Temp.: 30  $^{\circ}\text{C}$   
Detection: Suppressed Conductivity, Dionex ACES™  
Anion Capillary Electrolytic Suppressor  
Sample Prep.: 1:10 dilution

Peaks:	1. Fluoride	2.0 mg/L
	2. Chloride	10.0
	3. Nitrite	10.0
	4. Nitrate	10.0
	5. Bromide	10.0
	6. Sulfate	10.0
	7. Phosphate	20.0




29078

**FIGURE 5. Fast determination of inorganic anions in a municipal wastewater sample.**

Column: Dionex IonSwift MAX-200,  
0.25 mm × 250 mm  
Eluent Source: Dionex capillary EGC-KOH cartridge  
Gradient: KOH gradient from 2 mM to 50 mM  
Flow Rate: A: 10, B: 20, C: 25 µL/min  
Inj. Volume: 0.4 µL  
Column Temp.: 30 °C  
Detection: Suppressed Conductivity,  
Dionex ACES™ Anion Capillary  
Electrolytic Suppressor,  
Sample Prep.: 1:10 dilution

Peaks: 1. Chloride 25.5 mg/L  
2. Nitrate 73.3  
3. Sulfate 10.0




29079

**FIGURE 6. Fast separation of anions in a sports beverage by increasing flow rates.**

Column: Dionex IonSwift MAX-200,  
0.25 mm × 250 mm  
Eluent Source: Dionex capillary EGC-KOH cartridge  
Gradient: KOH gradient from 2 mM to 50 mM  
Flow Rate: A: 10, B: 20, C: 25 µL/min  
Inj. Volume: 0.4 µL  
Column Temp.: 30 °C  
Detection: Suppressed Conductivity, Dionex ACES™  
Anion Capillary Electrolytic Suppressor  
Sample Prep.: 1:10 dilution

Peaks: 1. Chloride 41.3 mg/L  
2. Phosphate 34.5  
3. Citrate 331



29080


## Oxyhalides Separation on the Dionex IonSwift MAX-200 Capillary Column

The Dionex IonSwift MAX-200 capillary column was selected for this application because the monolith backbone is ideal for fast separations using a hydroxide gradient, and the column chemistry is optimized for oxyhalide separations such as bromate, chlorite, and chlorate. The separation of bromate in municipal drinking water is shown Figure 7. Typically, municipal wastewater samples contain either chlorite and chlorate from chlorination or bromate from ozonation as a result from the disinfection processes. The samples rarely contain both types. However, this wastewater sample selected for bromate-spiking experiments already contained chlorite from the disinfectant treatment. In Figure 8, the chromatogram demonstrates that chlorite and chlorate can be separated from other inorganic anions using the Dionex IonSwift MAX-200 capillary column.

**FIGURE 7. Determination of bromate in a municipal drinking water.**

Column: Dionex IonSwift MAX-200, 0.25 × 250 mm  
 Eluent Source: Dionex EGC-KOH capillary  
 Gradient: 2 mM KOH from -10 to 0.1 min, 2–10mM from 0.1 to 10 min, 10–50 mM from 10 to 15 min, 50 mM from 15 to 20 min  
 Column Temp.: 30 °C  
 Flow Rate: 10 µL/min  
 Inj. Volume: 3 µL  
 Detection: Suppressed conductivity, Dionex ACES™ Anion Capillary Electrolytic Suppressor, recycle  
 Sample Prep.: 10 µg/L of bromate was spiked into a municipal drinking water  
 Peaks:


1. Fluoride	0.84 mg/L
2. Chlorite	—
3. Bromate	0.01
4. Chloride	4.05
5. Nitrate	0.21
6. Sulfate	2.14
7. Phosphate	0.07



**FIGURE 8. Determination of chlorite and chlorate in anion standard solution.**

Column: Dionex IonSwift MAX-200, 0.25 × 250 mm  
 Eluent Source: Dionex EGC-KOH capillary  
 Gradient: 2 mM KOH from -10 to 0.1 min, 2–10 mM from 0.1 to 10 min, 10–50 mM from 10 to 15 min, 50 mM from 15 to 20 min  
 Column Temp.: 30 °C  
 Flow Rate: 10 µL/min  
 Inj. Volume: 3 µL  
 Detection: Suppressed conductivity, Dionex ACES™ Anion Capillary Electrolytic Suppressor, recycle  
 Sample Prep.: 10 µg/L of bromate was spiked into a municipal drinking water  
 Peaks:

1. Fluoride	0.1 mg/L
2. Chlorite	0.1
3. Chloride	0.5
4. Nitrite	0.5
5. Chlorate	0.1
6. Bromide	0.5
7. Nitrate	0.5
8. Sulfate	0.5
9. Phosphate	1.0



29082

## Linearity and Detection Limit Determination

Calibration linearity for bromate was investigated by injecting standard solutions at four different concentrations (Table 1). Good linearity was observed from 3 to 30 µg/L for bromate. The method detection limit (MDL) was estimated to be 1 µg/L by examining the signal-to-noise ratio (3x S/N).

**Table 1. Linearity and MDL of Bromate Using the Dionex IonSwift MAX-200 Capillary Column at 10 µL/min Flow Rate Conditions**

Concentration of Standards (µg/L)	Correlation Coefficient (r <sup>2</sup> )	Calculated MDL (µg/L)
3	> 0.99	1.0
10		
20		
30		

## Recovery of Bromate in Municipal Drinking Water

To demonstrate the accuracy of the method, 10 µg/L of bromate was spiked into a municipal drinking water. An average of 87% recovery of bromate with 7% RSD was achieved. Due to co-elution of bromate and chlorite, the integration of bromate was affected resulting in low recovery. However, bromate and chlorite will not be present at the same time in a real municipal drinking water sample because the water is either treated by chlorination, generating chlorite and chlorate, or by ozonation, generating bromate.

## Conclusion

The Dionex IonSwift MAX-200 capillary column is a high-speed, high-resolution monolith column designed for the separation of oxyhalides and inorganic anions.<sup>1</sup> This column provides fast separation and high throughput. The high-pressure compatibility of the capillary IC format allows higher system backpressures further benefiting the IC capillary system to run fast, high-resolution ion analysis. Two-dimensional IC with suppressed conductivity or mass spectroscopy detection can be an alternative approach to determine oxyhalides if lower LOD is required.<sup>2,3</sup>

## References

1. Thermo Fisher Scientific. *Determination of Trace Concentrations of Chlorite, Bromate, and Chlorate in Bottled Natural Mineral Waters*, Dionex Application Note 184 (LPN 1890), Sunnyvale, CA, 2007.
2. Thermo Fisher Scientific. *Determination of Sub-µg/L Bromate in Municipal and Natural Mineral Waters Using Preconcentration with Two-Dimensional Ion Chromatography and Suppressed Conductivity Detection*, Dionex Application Note 187 (LPN 1943-02), Sunnyvale, CA, 2007.
3. Liu, Y.; Baretto, V.; Slingsby, R.; Pohl, C. *Determination of Trace-Level Bromate and Perchlorate in Environmental Samples Using Two-Dimensional Capillary Ion Chromatography with Suppressed Conductivity and Mass Spectroscopy Detection*. Thermo Fisher Scientific, Dionex Products Poster (LPN 2243), Sunnyvale, CA, 2009.

©2011 Thermo Fisher Scientific, Inc. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others.

**Dionex Products:** 1228 Titan Way, PO Box 3603, Sunnyvale, CA 94088-3603, (408) 737-0700  
**North America:** U.S./Canada (847) 295-7500  
**South America:** Brazil (55) 11 3731 5140  
**Europe:** Austria (43) 616 51 25, Benelux (31) 20 683 9768 (32) 3 353 4294  
Denmark (45) 36 36 90 90, France (33) 1 39 30 01 10, Germany (49) 61125 991 0  
Ireland (353) 644 0064, Italy (39) 02 51 62 1267, Sweden (46) 8 473 3380,  
Switzerland (41) 62 205 9966, United Kingdom (44) 1276 691722  
**Asia Pacific:** Australia (61) 2 9420 5233, China (852) 2428 3282, India (91) 22 2764 2735,  
Japan (81) 6885 1213, Korea (82) 2 2653 2580, Singapore (65) 6289 1190,  
Taiwan (886) 2 875 6655