

# Automated Solid-Phase Extraction of Wastewaters and Surface Waters Using the AutoTrace<sup>®</sup> SPE Instrument for EPA Method 608

## **INTRODUCTION**

EPA method 608 uses liquid-liquid extraction by Separatory funnel combined with gas chromatography and electron capture detector for the analysis of organochlorine pesticides and PCBs in municipal and industrial waste water.

Unlike the drinking water methods, however, the extraction of wastewaters may contain particulates or suspensions that may block or clog SPE materials. The AutoTrace SPE instrument has features that allow extraction of samples even with high amounts of suspended solids.

The loading pumps are able to load samples at a given rate even as the bed begins to block. If the bed does get completely blocked, the sample is still in the original container from where it can be retrieved for further processing without loss. The solvent delivery syringe is able to deliver solvent under sufficient positive pressure to overcome severe blockages.

## **AUTOTRACE EXTRACTION PROCEDURES**

2.5 mL of methanol is added to 500 mL of sample. Various river water samples (500 mL) were spiked 100  $\mu$ L of surrogate with 5ng/ $\mu$ L each of Tetrachlorometaxylene (TCMX) and Decachlorobiphenyl (DCB). Analytes were spiked with 100  $\mu$ L at 2.5 ng/ $\mu$ L.

### **SPE Cartridge Used:**

C18 cartridge 500 mg in a 6 mL cartridge with a polyurethane plug prefilter, part number 221-0050-CD.



*AutoTrace Instrument for Automated SPE*

## **AUTOTRACE EXTRACTION PROCEDURE**

### **Process 6 Samples using the following procedure:**

- Step 1: Wash syringe with 2 mL of CH<sub>3</sub>OH
- Step 2: Rinse column with 5 mL of EtOAc into solvent waste
- Step 3: Rinse column with 5 mL of DCM into solvent waste
- Step 4: Condition column with 10 mL of CH<sub>3</sub>OH into solvent waste
- Step 5: Condition column with 10 mL of water into aqueous waste
- Step 6: Pause and alert operator, resume when Continue is pressed
- Step 7: Load 550 mL of sample onto column
- Step 8: Dry column with gas for 10 minutes

| <b>Solvents</b>               | <b>Flow Rate</b>                                  |
|-------------------------------|---|
| Water (reagent grade)         | Cond. Flow: 40 mL/min                             |
| CH <sub>3</sub> OH (methanol) | Load Flow: 3 mL/min                               |
| EtOAc (ethyl acetate)         | Rinse Flow: 40 mL/min                             |
| DCM (dichloromethane)         | Cond. Air Push: 15 mL/min                         |
| Unused (load with methanol)   | Elute Flow: 20 mL/min<br>Elute Air Push: 5 mL/min |

Step 6 is added so the conditioning steps can be started before final sample preparations are completed. This alleviates any concern that the instrument will load the sample before the operator is ready.

## **SAMPLE ELUTE PROCEDURE**

### **Process 6 Samples using the following procedure:**

- Step 1: Collect 5 mL fraction into sample tube using EtOAc
- Step 2: Manually Rinse sample container with 7 mL to collect
- Step 3: Manually Rinse sample container with 10 mL to collect
- Step 4: Collect 2 mL fraction into sample tube using DCM

| <b>Solvents</b>               | <b>Flow Rate</b>                                  |
|-------------------------------|---|
| Water (reagent grade)         | Cond. Flow: 40 mL/min                             |
| CH <sub>3</sub> OH (methanol) | Load Flow: 3 mL/min                               |
| EtOAc (ethyl acetate)         | Rinse Flow: 40 mL/min                             |
| DCM (dichloromethane)         | Cond Air Push: 15 mL/min                          |
| Unused (load with methanol)   | Elute Flow: 20 mL/min<br>Elute Air Push: 5 mL/min |

*Note change in Load Flow rate.*

### **Manual Rinse Details**

At Step 2, add 5 mL of ethyl acetate to sample container; at Step 3 add 5 mL of dichloromethane to sample container.

All channels that succeed in delivering solvent from their sample containers are now finished. The channel's tension is released and the extracts removed, dried, and concentrated to an appropriate final volume as dictated by method 608 or Detection Limit requirements.

Wastewaters and surface water often contain high levels of suspended solids. The selection of ISOLUTE® Depth Filters can trap particulate above the reservoir. The AutoTrace system can process up to six samples. If any sample plugs the SPE Cartridge, simply release the tension on that channel to prevent further processing and continue to process the remaining samples during sample loading.

For any Channel that does not push the contents through the SPE cartridge, simply load the Blocked Shocker program and run the program on the affected channel.

---

### **BLOCKED SHOCKER PROGRAM:**

#### **Process 6 Samples using the following procedure:**

- Step 1: Collect 3 mL fraction into sample tube using EtOAc
- Step 2: Manually Rinse sample container with 15 mL to collect
- Step 3: Collect 2 mL fraction into sample tube using DCM

Do not add more solvent at Step 3. The purpose of this program is to empty the sample container of the rinse solvent.

If loading pump is still unable to deliver solvent through SPE cartridge, the Blocked Shocker program can be run twice more. The solvent forced through by the piston will suffice to extract the analytes from the SPE cartridge.

Once all channels that successfully loaded have been eluted, resume work on samples that blocked by starting at the Sample Load Program with a fresh SPE cartridge on a different channel, loading the remainder of the sample. Elute both cartridges as described above and combine the extracts into one.

### **DISCUSSION**

The solvent delivery piston is more able to force liquid through the SPE than is the loading pump. Sometimes, the resistance is caused by residual water surface tension in the SPE bed and forcing some solvent through with the syringe will enable the resistance to lessen for the loading pump. If the sample contains sufficient solids to plug an SPE cartridge, the likelihood that analytes have adsorbed to the wall of the container and not been transferred to the cartridge is reduced, thus the solvent delivered by Blocked Shocker can be assumed to contain all of the analytes present in the sample.

Table 1 shows recoveries from a run (n=4) where all samples eluted normally. Table 2 shows recoveries from another run (n=4) that experienced various problems. Replicate S4788.2 resisted elution by the loading pumps and was eluted with 3 repetitions of Blocked Shocker. Replicates S4788.3 and S4788.5 needed Blocked Shocker once and then succeeded in eluting using the loading pumps. Replicate S4788.6 blocked before loading all of the sample and was restarted on a second SPE cartridge, then both elutions were combined.

#### **Analysis of Example Spike Recoveries**

Samples were analyzed on a gas chromatograph equipped with electron capture detectors. The following results were obtained using an Agilent® 6890 with dual 0.32 mm ID, 30 m long capillary columns, Rtx®-CLPesticides and Rtx-CLPesticidesII from Restek®.

**Table 1. Iowa River Water, Estimated Total Suspended Solids, 70 mg/L**

|                           | <b>S4785.2</b> | <b>S4785.3</b> | <b>S4785.5</b> | <b>S4785.6</b> | <b>Mean</b> | <b>Std. Dev.</b> |
|---------------------------|----------------|----------------|----------------|----------------|-------------|------------------|
| <b>SS:TCX</b>             | 72.7%          | 68.0%          | 68.3%          | 76.1%          | 71.3%       | 3.9%             |
| <b>SS:DCBP</b>            | 44.5%          | 70.2%          | 58.9%          | 48.6%          | 55.5%       | 11.5%            |
| <b>alpha-BHC</b>          | 98.7%          | 96.0%          | 98.3%          | 99.9%          | 98.2%       | 1.6%             |
| <b>gamma-BHC</b>          | 108.5%         | 119.2%         | 114.4%         | 110.9%         | 113.3%      | 4.7%             |
| <b>beta-BHC</b>           | 107.4%         | 158.0%         | 122.4%         | 112.2%         | 125.0%      | 22.9%            |
| <b>delta-BHC</b>          | 108.1%         | 150.0%         | 122.0%         | 112.7%         | 123.2%      | 18.8%            |
| <b>Heptachlor</b>         | 86.5%          | 93.0%          | 91.8%          | 89.3%          | 90.1%       | 2.9%             |
| <b>Aldrin</b>             | 63.2%          | 74.8%          | 71.8%          | 69.8%          | 69.9%       | 4.9%             |
| <b>Hepta. epoxide</b>     | 94.4%          | 124.6%         | 97.8%          | 99.7%          | 104.1%      | 13.8%            |
| <b>DDE</b>                | 67.3%          | 96.0%          | 78.9%          | 68.3%          | 77.7%       | 13.3%            |
| <b>Endosulfan I</b>       | 102.8%         | 145.3%         | 118.8%         | 110.3%         | 119.3%      | 18.5%            |
| <b>Dieldrin</b>           | 105.5%         | 152.6%         | 122.8%         | 113.3%         | 123.6%      | 20.6%            |
| <b>Endrin</b>             | 104.6%         | 151.8%         | 119.3%         | 108.6%         | 121.1%      | 21.4%            |
| <b>DDD</b>                | 87.6%          | 141.1%         | 105.3%         | 93.1%          | 106.8%      | 24.0%            |
| <b>Endosulfan II</b>      | 99.6%          | 158.1%         | 116.4%         | 105.0%         | 119.8%      | 26.5%            |
| <b>DDT</b>                | 80.5%          | 90.1%          | 65.6%          | 52.1%          | 72.1%       | 16.7%            |
| <b>Methoxychlor</b>       | 129.6%         | 177.7%         | 108.8%         | 86.9%          | 125.8%      | 38.8%            |
| <b>Endosulfan sulfate</b> | 122.0%         | 198.5%         | 142.6%         | 128.0%         | 147.8%      | 34.9%            |

**Table 2. Cedar River Water, Estimated Total Suspended Solids, 64 mg/L**

|                           | <b>S4788.2</b> | <b>S4788.3</b> | <b>S4788.5</b> | <b>S4788.6</b> | <b>Mean</b> | <b>Std. Dev.</b> |
|---------------------------|----------------|----------------|----------------|----------------|-------------|------------------|
| <b>SS:TCX</b>             | 54.9%          | 46.8%          | 49.1%          | 54.9%          | 51.4%       | 4.1%             |
| <b>SS:DCBP</b>            | 23.6%          | 19.0%          | 25.6%          | 33.2%          | 25.4%       | 5.9%             |
| <b>alpha-BHC</b>          | 64.6%          | 48.4%          | 58.8%          | 91.2%          | 65.7%       | 18.3%            |
| <b>gamma-BHC</b>          | 72.1%          | 53.7%          | 64.8%          | 97.6%          | 72.0%       | 18.6%            |
| <b>beta-BHC</b>           | 68.5%          | 54.1%          | 63.9%          | 94.5%          | 70.2%       | 17.2%            |
| <b>delta-BHC</b>          | 69.0%          | 53.0%          | 64.3%          | 97.2%          | 70.9%       | 18.8%            |
| <b>Heptachlor</b>         | 49.0%          | 36.7%          | 38.6%          | 59.0%          | 45.8%       | 10.3%            |
| <b>Aldrin</b>             | 40.2%          | 30.7%          | 35.6%          | 50.9%          | 39.4%       | 8.6%             |
| <b>Hepta.epoxide</b>      | 62.5%          | 48.6%          | 58.3%          | 81.8%          | 62.8%       | 13.9%            |
| <b>DDE</b>                | 37.3%          | 28.8%          | 34.2%          | 49.5%          | 37.5%       | 8.8%             |
| <b>Endosulfan I</b>       | 70.7%          | 54.9%          | 65.3%          | 87.7%          | 69.7%       | 13.7%            |
| <b>Dieldrin</b>           | 70.1%          | 55.0%          | 65.9%          | 90.4%          | 70.4%       | 14.8%            |
| <b>Endrin</b>             | 63.1%          | 48.6%          | 55.3%          | 77.5%          | 61.1%       | 12.4%            |
| <b>DDD</b>                | 50.8%          | 38.7%          | 46.2%          | 62.5%          | 49.6%       | 10.0%            |
| <b>Endosulfan II</b>      | 68.6%          | 52.7%          | 64.1%          | 79.0%          | 66.1%       | 10.9%            |
| <b>DDT</b>                | 23.5%          | 17.6%          | 19.4%          | 26.5%          | 21.8%       | 4.0%             |
| <b>Methoxychlor</b>       | 47.7%          | 35.5%          | 32.7%          | 65.2%          | 45.3%       | 14.8%            |
| <b>Endosulfan sulfate</b> | 83.6%          | 63.8%          | 76.8%          | 101.7%         | 81.5%       | 15.8%            |

**ACKNOWLEDGEMENTS**

Data and study submitted by: David Larabee-Zierath,  
PhD of the University Hygienic Laboratory, Iowa's  
Environmental and Public Health Laboratory, Iowa City, IA

AutoTrace is a registered trademark of Dionex.  
ISOLUTE is a registered trademark of Biotage.  
Agilent 6890 is a registered trademark of Agilent.  
Rtx and Restek are registered trademarks of the Restek Corporation.

**Passion. Power. Productivity.**



**Dionex Corporation**  
1228 Titan Way  
P.O. Box 3603  
Sunnyvale, CA  
94088-3603  
(408) 737-0700

**North America**  
U.S. (847) 295-7500  
Canada (905) 844-9650  
**South America**  
Brazil (55) 11 3731 5140

**Europe**  
Austria (43) 1 616 51 25 Benelux (31) 20 683 9768 (32) 3 353 4294  
Denmark (45) 36 36 90 90 France (33) 1 39 30 01 10 Germany (49) 6126 991 0  
Ireland (353) 1 644 0064 Italy (39) 02 51 62 1267 Sweden (46) 8 473 3380  
Switzerland (41) 62 205 9966 United Kingdom (44) 1276 691722

**Asia Pacific**  
Australia (61) 2 9420 5233 China (852) 2428 3282 India (91) 22 2764 2735  
Japan (81) 6 6885 1213 Korea (82) 2 2653 2580 Singapore (65) 6289 1190  
Taiwan (886) 2 8751 6655

[www.dionex.com](http://www.dionex.com)



LPN 2163 PDF 01/09  
©2009 Dionex Corporation