

# Automated Solid Phase Extraction (SPE) of Dioxins and Furans in Surface Water

## **INTRODUCTION**

The EPA Methods 8260 and 8270 are the basis of the sample preparation of polychlorinated dibenzo-p-dioxins, polychlorinated dibenzofurans and PCBs in surface and ground water. The clean up described in this Application Brief (AB) also recovers 1,4-Dioxan. The addition of a 1% v/v isopropanol solution to the sample bottle reduces losses due to adhesion to glass and tubing. The use of HCL 0.1 M to acidify the sample to pH 2 prevents the humic acids often present in the water from allowing dioxins to pass through the C18 cartridge.

## **INSTRUMENTATION USED FOR SAMPLE PREPARATION**

Dionex AutoTrace® instrument

## **SAMPLE PREPARATION**

- A 500 mL sample of river water is taken for the analysis and mixed with 5 mL of isopropanol to prevent the components from sticking to the glass bottle.
- The use of HCL 0.1 M is important to adjust the water to pH 2 to enhance recoveries
- A standard C18 end capped SPE silica column, 3 mL, 1g, can be used for the clean up. The elution solvent is ethanol/toluene (70:30) mL and six sample collection positions for the eluent.

## **AUTOTRACE METHOD AND CONCLUSION**

All solvent lines are purged and primed with solvent first. A maximum of five solvents can be used to run a wide range of methods. Independent lines separate the aqueous and organic waste solvents. The instrument sample rack has six positions that can take volumes from 10 mL to 2000 mL and six sample collection positions for the eluent. This method offers an automated rugged and reproducible solution for cleaning up the samples to concentrate and remove interferences.

<b>No.</b>	<b>Method: Estimated time 1h 30 min</b>
1	Process six samples using the following method steps:
2	Condition column with 5.0 mL of CH <sub>3</sub> OH into solvent waste.
3	Condition column with 5.0 mL of water into aqueous waste.
4	Load 500.0 mL of sample into column.
5	Rinse column with 3.0 mL of CH <sub>3</sub> OH/water 40:60 into solvent waste.
6	Wash syringe with 5.0 mL of 70%EtOH/toluene.
7	Dry column with gas for 10 min.
8	Collect 5.0 mL fraction into sample tube using 70%EtOH/toluene.
9	Wash syringe with 5.0 mL of CH <sub>3</sub> OH.
10	End

---

**PARAMETERS:****Flow Rates**

Cond Flow:	15.0 mL/min
Load Flow:	10.0 mL/min
Rinse Flow:	20.0 mL/min
Elute Flow:	5.0 mL/min
Cond Air Push:	15.0 mL/min
Rinse Air Push:	20.0 mL/min
Elute Air Push:	5.0 mL/min

**SPE Parameters**

Push Delay:	5 sec
Air Factor:	1.0
Autowash Vol.:	1.00 mL

**Instrument Parameters**

Max. Elution Vol.:	12.0 mL
Exhaust Fan On:	Yes
Beeper On:	Yes

AutoTrace is a registered trademark of Dionex Corporation.

Passion. Power. Productivity.

**Dionex Corporation**

1228 Titan Way  
P.O. Box 3603  
Sunnyvale, CA  
94088-3603  
(408) 737-0700

**North America**

U.S. (847) 295-7500

**South America**

Brazil (55) 11 3731 5140

**Europe**

Austria (43) 1 616 51 25 Benelux (31) 20 683 9768; (32) 3 353 4294  
Denmark (45) 36 36 90 90 France (33) 1 39 30 01 10 Germany (49) 6126 991 0  
Ireland (353) 1 644 0064 Italy (39) 02 51 62 1267 Sweden (46) 8 473 3380  
Switzerland (41) 62 205 9966 United Kingdom (44) 1276 691722

**Asia Pacific**

Australia (61) 2 9420 5233 China (852) 2428 3282 India (91) 22 2764 2735  
Japan (81) 6 6885 1213 Korea (82) 2 2653 2580 Singapore (65) 6289 1190  
Taiwan (886) 2 8751 6655

[www.dionex.com](http://www.dionex.com)



LPN 2214 PDF 02/09  
©2009 Dionex Corporation