

# Automated Solid Phase Extraction (SPE) of Total Petroleum Hydrocarbons Using Dionex AutoTrace® Instrument

## **INTRODUCTION**

This method was developed for the extraction of total petroleum hydrocarbons (or silica gel treated-hexane extractable material [SGT-HEM] as described in EPA Method 1664) from water samples. The analytes are determined by gravimetric analysis.

## **Extraction Procedure**

**\*ISOLUTE® SPE Column:** TPH 1 g / 6 ml

**Pre-treatment:** Add 10 mL of methanol to 1 L of sample. Acidify the sample with 6 M HCl to a pH ~2.

**Solvation:** Program the AutoTrace instrument to WASH SYRINGE with 5 mL of methanol, followed by a CONDITION COLUMN with 10 mL of methanol.

**Equilibration:** Program the AutoTrace instrument to WASH SYRINGE with 10 mL of distilled deionised water pH~2, followed by a CONDITION COLUMN with 10 mL of distilled deionised water pH~2.

**Sample application:** Program AutoTrace instrument to LOAD 1100 mL SAMPLE. Following the sample load step, include a PAUSE AND ALERT step.

**Interference elution:** During pause step, remove sample line from sample bottle. Manually rinse sample bottle with 10 mL of acetone, swirl well to cover sides of bottle. Dilute acetone in bottle with 40 mL of distilled deionised water, pH~2. Return sample line to bottle, and press CONT. Program the AutoTrace to LOAD 60 mL SAMPLE. Following this step, include a PAUSE AND ALERT step. Repeat this step if necessary until the rinse solvent appears clear. Press CONT. Dry cartridge assembly for 30 min.

**Analyte elution:** Rinse the AutoTrace syringe with 5 mL of hexane. ELUTE TO SOAK with 4 mL of hexane. Follow this step with a TIMED PAUSE for 2 min. ELUTE TO COLLECT with 4 mL of hexane into tared collection tubes. Concentrate solvent to near dryness in a gentle stream of nitrogen. Collection tubes may be placed in a heating block held at 35 °C to expediate evaporation. Once solvent is almost gone, weigh tubes in one minute intervals until weight loss is less than 1 mg. Total petroleum hydrocarbons is the weight of residue.

## **Setup parameters on AutoTrace should be as follows:**

### **Flow Rates**

Cond Flow: 40.0 ml/min

Cond Air Push: 40.0 ml/min

Load Flow: 30.0 ml/min

Rinse Air Push: 40.0 ml/min

Rinse Flow: 40.0 ml/min

Elute Air Push: 40.0 ml/min

Elute Flow: 20.0 ml/min

### **SPE Parameters**

Push Delay: 5 sec

Air Factor: 1.0

Autowash Vol.: 1.00 ml

\* ISOLUTE column part numbers represent the product configuration of choice for use with the AutoTrace instrument.

**THE AUTOTRACE INSTRUMENT PROCEDURE SHOULD BE WRITTEN AS FOLLOWS:**

Step	Method:
1	Process 6 Samples using the following method steps:
2	Wash syringe with 5 mL of CH <sub>3</sub> OH.
3	Condition column with 10 mL of CH <sub>3</sub> OH.
4	Wash syringe with 10 mL of water pH=2.
5	Condition column with 10 mL of water to aqueous waste.
6	Load 1100.0 mL of sample on column.
7	Pause and Alert operator, resume when CONTInue is pressed.
8	Load 60 mL of sample onto column.
9	Pause and Alert operator, resume when CONTInue is pressed.
10	Dry column with gas for 30 min.
11	Wash syringe with 5 mL of hexane.
12	Soak and collect 4 mL fraction using hexane.
13	Pause for 2 min.
14	Collect 3 mL fraction into sample tube using hexane.
15	Clean each sample path with 10 mL into solvent waste.
16	Clean each sample path with 50 mL into aqueous waste.
17	Wash syringe with 10 mL of water.
18	End

**Structural Various:**

Non-polar hydrocarbons.

**Structural Considerations:**

This method is suitable for the non-polar petroleum hydrocarbons (SGT-HEM).

**Matrix Considerations:**

The matrix is polar, and the analytes are extracted by a non-polar retention mechanism.

**Analytical Method:**

Gravimetric analysis using an analytical balance having precision to 0.1 mg.

**General Comments:**

- This method describes an automated procedure for the determination of silica gel treated-hexane extractable material (SGT-HEM) from an aqueous sample. This column can also be used for total oil and grease type extractions (EPA 1664 HEM) if the concentration of polar components is low. In this situation, the hexane elution step should be followed by an elution (into a second tared vial) using 2 x 4 mL THF/hexane (1:1, v/v). Combination of these two fractions will give a determination of total oil and grease, as described in EPA Method 1664. However, for higher concentrations of polar components, refer to application note IST1018 (TPH/Total Oil and Grease determination on the AutoTrace instrument) which uses a column specifically optimized for EPA 1664 extractions on the AutoTrace instrument.
- Due to the nature of the analytes, the bottle washing steps after sample loading are very important, as analytes do stick to the walls of the sample bottle. For this reason, sample splitting is not recommended for TPH or oil and grease samples.

**CONCLUSION**

The AutoTrace instrument is specifically designed for extraction of large volume water samples. The extraction of silica gel treated hexane extractable material (SGT-HEM) is efficient using Isolute SPE column TPH. The positive flow feature of the AutoTrace instrument ensures that each extraction is uniform and optimizes EPA Method 1664.

AutoTrace is a registered trademark of Dionex Corporation.

Passion. Power. Productivity.



**Dionex Corporation**

1228 Titan Way  
P.O. Box 3603  
Sunnyvale, CA  
94088-3603  
(408) 737-0700

**North America**

U.S. (847) 295-7500

**South America**

Brazil (55) 11 3731 5140

**Europe**

Austria (43) 1 616 51 25 Benelux (31) 20 683 9768; (32) 3 353 4294  
Denmark (45) 36 36 90 90 France (33) 1 39 30 01 10 Germany (49) 6126 991 0  
Ireland (353) 1 644 0064 Italy (39) 02 51 62 1267 Sweden (46) 8 473 3380  
Switzerland (41) 62 205 9966 United Kingdom (44) 1276 691722

**Asia Pacific**

Australia (61) 2 9420 5233 China (852) 2428 3282 India (91) 22 2764 2735  
Japan (81) 6 6885 1213 Korea (82) 2 2653 2580 Singapore (65) 6289 1190  
Taiwan (886) 2 8751 6655

[www.dionex.com](http://www.dionex.com)



LPN 2226 PDF 02/09  
©2009 Dionex Corporation