

Rapid HPLC Separation of Multiclass Antibiotics in Food and Water

INTRODUCTION

The use of antibiotics on livestock, aquaculture, and bee husbandry helps maintain health and provides other benefits such as improved disease resistance, increased production, and in some cases, reduction in foodborne pathogens. However, antibiotic residues in foods can cause undesirable side effects such as idiosyncratic aplastic anemia, production of antibiotic-resistant bacteria, and the reduction of indigenous microbiota found in the human digestive tract. In addition, the wastes of these animals can lead to antibiotics being present in our water. Organizations responsible for food safety perform assays for the presence of allowed and banned antibiotics as part of their surveillance activities. This application brief demonstrates a reproducible Rapid Separation Liquid Chromatography (RSLC) method for the separation of polyketide (tetracycline and oxytetracycline), macrolide (tylosin), phenicol (chloramphenicol), nitrofurantoin, and sulfonamide (sulfathiazole) classes of antibiotics with an MS-compatible mobile phase (Figure 1, Table 1).

Now sold under the
Thermo Scientific brand





Figure 1. Separation of multiple classes of antibiotics using the Dionex RSLC system and UV-diode array detection at 267 nm.

CONDITIONS

System: UltiMate® 3000 RSLC consisting of HPG-3200 Binary pump, SRD-3200 solvent rack with degasser channels, WPS-3000TRS thermostatted in-line split-loop autosampler, TCC-3000RS thermostatted column compartment, and DAD-3000RS diode array detector with semi-micro (2.5 µL, 7 mm) SS flow cell

Column: Acclaim® RSLC Polar Advantage II (2.1 x 150 mm packed with 2.2 µm particles, 120 Å average pore size), P/N 071401

Mobile Phases: A = 0.5% formic acid in water (pH 2.4)

B = 0.5% formic acid in acetonitrile

Gradient:	Time (min)	%A	%B
	0.00	90.0	10.0
	2.00	74.0	26.0
	4.00	48.0	52.0
	4.10	10.0	90.0
	5.90	10.0	90.0
	6.00	90.0	10.0

Maximum pressure – 685 bar (9900 psi)

Flow Rate: 0.80 mL/min

Inj. Volume: 10.0 µL

Temperatures: 30 °C (column)

5 °C (autosampler)

Detection: 267, 278, 287, 355, 367 nm (based on UV spectra of individual antibiotics)

Bandwidth – 11 nm each

Data collection rate – 25 Hz

Response time – 0.20 s

Table 1. Retention Time Precision for the Separation of Multiclass Antibiotics^{a, b}

Peak	Antibiotic	RT (min)	RSD
1	Sulfathiazole	1.725	0.13
2	Oxytetracycline	1.925	0.11
3	Tetracycline	2.109	0.12
4	Nitrofurantoin	2.305	0.12
5a	Tylosin & Related Products	3.492	0.12
5b		3.673	0.11
5c		3.789	0.11
6	Chloramphenicol	3.610	0.12

^a Acclaim RSLC Polar Advantage II column

^b n = 10 replicates

Acclaim and UltiMate are registered trademarks of Dionex Corporation.

Passion. Power. Productivity.



Dionex Corporation

1228 Titan Way
P.O. Box 3603
Sunnyvale, CA
94088-3603
(408) 737-0700

North America

U.S./Canada (847) 295-7500

South America

Brazil (55) 11 3731 5140

Europe

Austria (43) 1 616 51 25 Benelux (31) 20 683 9768 (32) 3 353 4294
Denmark (45) 36 36 90 90 France (33) 1 39 30 01 10 Germany (49) 6126 991 0
Ireland (353) 1 644 0064 Italy (39) 02 51 62 1267 Sweden (46) 8 473 3380
Switzerland (41) 62 205 9966 United Kingdom (44) 1276 691722

Asia Pacific

Australia (61) 2 9420 5233 China (852) 2428 3282 India (91) 22 2764 2735
Japan (81) 6 6885 1213 Korea (82) 2 2653 2580 Singapore (65) 6289 1190
Taiwan (886) 2 8751 6655

www.dionex.com



LPN 2412 PDF 01/10
©2010 Dionex Corporation