

New Developments in Capillary Ion Chromatography Systems with Electrochemical Detection and Their Applications

Now sold under the
Thermo Scientific brand

Thermo
SCIENTIFIC

Jun Cheng, Petr Jandik, Chris Pohl, and Jason Wood, Dionex Corporation, Sunnyvale, CA, U.S.A

ABSTRACT

There has been increasing interest in the development of capillary ion chromatography (IC) systems and methods for the determination of electrochemically active species. The practice of IC in the capillary format (i.e., using small-bore columns with internal diameters of < 1 mm) has a number of potential advantages including small sampling volumes, potential of improved mass detection limits, and lower eluent consumption. Here, the use of such a system in carbohydrates analysis is presented. Low sample volumes, in particular, offer improved compatibility with Bio IC applications, where biologically relevant samples are limited. Additionally, the operation of capillary IC at low flow rates reduces waste generated in the laboratory setting. The practice of IC in the capillary format offers possibilities of new selectivity for difficult applications using new columns packed with previously difficult-to-prepare stationary phases.

Here, ongoing efforts in the development of capillary IC systems with on-line electrolytic eluent generation (EG) and electrochemical detection for determination of target analytes are reported. A discussion is presented on the optimization of capillary consumables and the advantages of using capillary IC systems in the determination of target analytes. Furthermore, it is also demonstrated that with *Always on, Always Ready*, a system at capillary chemistries and flow rates allows undiminished sensitivity and selectivity with minimum involvement by the IC user.

INTRODUCTION

Ion chromatography is a prominent technique to determine inorganic and organic ions in all types of liquid matrices. Monitoring ions in drinking, surface, waste, and process waters is routinely performed using IC with either laboratory or process IC systems. Common applications include testing of disinfectant byproducts such as bromate in tap or bottled water, iodide and iodate in sea water, amines in wastewater, and trace-level anions and cations in ultrapure water for the semiconductor and power industry. It is important to determine the type and composition of carbohydrates in food, beverage, and dairy samples. Moreover, monosaccharide composition analysis is a powerful tool in the characterization of glycoproteins. The results are useful in glycoprotein drug development. Most of these applications are performed using column chemistries in 3 or 4 mm i.d. format and flow rates of 0.50–1 mL/min.

Here, the authors present the next evolution in IC using capillary chemistries. Capillary chemistries use 0.4 mm i.d. columns at flow rates of 10 μ L/min. This allows an IC system to be operated continuously for up to 12 months with only 5.2 L of water. Furthermore, less time is required for eluent preparation and equilibration or calibration because the system can be always ready for the next sample. Due to the low eluent consumption, the waste generated is also reduced, which makes capillary IC beneficial when toxic or radioactive samples are analyzed. Also, the injection volume is reduced to 400 nL typically while maintaining minimum detection limits. This makes capillary IC suitable for applications that are sample-limited, such as metabolomics, corrosion monitoring of individual corrosion spots, and process monitoring of small-scale reactors.

The Dionex ICS-5000 system is the first IC system that supports analytical scales (4 and 2 mm) and capillary chemistries (Figure 1). The dual system allows users to operate one channel with 2 or 4 mm column chemistries and explore capillary chemistries on the other channel. The ICS-5000 system supports manually prepared eluents as well as Reagent-Free™ IC technologies on analytical and capillary scales, and eluent regeneration for 4 mm column chemistries (Figure 2).



Figure 1. ICS-5000 system with IC Cube™ module.

The specifications of the Dionex ICS-5000 system are listed below:

Features	ICS-5000 System
Detector Compartment Temperature zones	Dual Zone Lower zone has one temperature; both cubes have independent temperature settings
Suppressors Capacity Void volumes	Anion Capillary Electrolytic Suppressor (ACES) 2 µeq/min < 1.0 µL
Pump Technology Pump flow rates Typical flow rates Maximum pressure Pressure ripple Flow accuracy	Dual-Piston (in series) 0.001– 3.000 mL/min 5–20 µL 6000 psi (capillary flow rates) < 0.25% at 10 µL/min < 0.1%
Detection Technology Cell volume (at working electrode) Reference electrode Noise specification	Electrochemical Detector < 130 nL Palladium hydrogen < 10 pA or 20 pC
Eluent Generation Concentration range Lifetime	RFIC-EG™ 0.1–200 mM (capillary) 18 months continuous operation (75 mM at 10 µL/min)




Figure 2. Block diagram of a capillary reagent-free ion chromatography (RFIC™) system with electrochemical detection.




Figure 3. PdH reference electrode.

SYSTEM INTEGRATION WITH THE IC CUBE

In liquid chromatography, often smaller i.d. columns (especially nano and capillary columns) can be difficult to use. Operators need to be careful in making fluidic connections and ensuring that they are dead-volume free. In addition, all components need to be volume optimized for low flow rates of 10 µL/min. In capillary scale, a dead volume of several nanoliters can broaden peaks and dramatically affect the chromatographic resolution. To overcome these issues, the ICS-5000 system has all capillary consumables in a separate compartment called the IC Cube™ module. The IC Cube module houses the capillary column and guard column in a cartridge, the injection valve with an internal loop of 400 nL, the capillary suppressor module, an optional carbonate removal device, and the eluent generator degasser module. The usage of cartridges simplifies the plumbing by reducing the number of fluidic connections by 50% compared to the same analytical configuration.



The IC Cube

RESULTS

The newly designed, capillary-based Dionex ICS-5000 system was used to detect monosaccharides in various matrices, including orange and cranberry-raspberry juice (Figures 4–8). The system was initially qualified with a mix of six monosaccharides and limits of detection (LOD) were determined (see Table 1).

After completing the initial runs on six monosaccharides, alditol, and disaccharides, the ICS-5000 system was used to determine glucose, fructose, and sucrose in two fruit juices to evaluate the system with actual samples (matrices).

2 New Developments in Capillary Ion Chromatography Systems with Electrochemical Detection and Their Applications

Table 1. Analytical Performance: Six Monosaccharides

	LOD (μM)	Linear Range (μM)	Correlation Coefficient
Fuc	0.024	0.024–50	0.9958
GalN	0.018	0.018–25	0.9998
GlcN	0.029	0.029–25	0.9949
Gal	0.054	0.054–25	0.9964
Glc	0.056	0.056–50	0.9986
Man	0.068	0.068–50	0.9956

Note: Limit of detection (LOD) is three times the noise; injection volume is 400 nL.




Figure 4. Separation of six monosaccharides.




Figure 5. Two-day detection stability of six monosaccharides.




Figure 6. Separation of mono- and disaccharides.




Figure 7. Separation of alditol and disaccharides.




Figure 8. Analysis of fruit juice samples.

SUMMARY

Here, the separation and detection of alditols, mono-, and disaccharides using capillary IC and electrochemical detection is demonstrated. The ICS-5000 system provides excellent limits of detection, reproducibility, and linear range. The long term stability of the system is excellent, with < 3% RSDs for area reproducibility. The advantages of the capillary system include:

- Minimal eluent consumption (< 15 mL/day)
- Improved detection limits for monosaccharides (< 10–70 nM)
- Small sample size (400 nL injection volumes)
- Capillary consumables integrated into the IC Cube for ease of access and use

CarboPac is a registered trademark and IC Cube, Reagent-Free, RFIC, and RFIC-EG are trademarks of Dionex Corporation.

Passion. Power. Productivity.



Dionex Corporation

1228 Titan Way
P.O. Box 3603
Sunnyvale, CA
94088-3603
(408) 737-0700

North America

U.S./Canada (847) 295-7500

South America

Brazil (55) 11 3731 5140

Europe

Austria (43) 1 616 51 25 Benelux (31) 20 683 9768 (32) 3 353 4294
Denmark (45) 36 36 90 90 France (33) 1 39 30 01 10 Germany (49) 6126 991 0
Ireland (353) 1 644 0064 Italy (39) 02 51 62 1267 Sweden (46) 8 473 3380
Switzerland (41) 62 205 9966 United Kingdom (44) 1276 691722

Asia Pacific

Australia (61) 2 9420 5233 China (852) 2428 3282 India (91) 22 2764 2735
Japan (81) 6 6885 1213 Korea (82) 2 2653 2580 Singapore (65) 6289 1190
Taiwan (886) 2 8751 6655

www.dionex.com



LPN 2444-01 2/10
©2010 Dionex Corporation